



# ETHOS X APPLICATION REPORTS

Microwave Green Extraction of Natural Products



**MILESTONE**  
HELPING  
CHEMISTS


# Index

- Introduction
- Important parameters in Microwave-Assisted Extraction
- Microwave Extraction Techniques
- How does Milestone conceive these extraction techniques?
- EX00 - BASIL FRAGRANCES
- EX01 - CITRUS PEEL FLAVORS
- EX02 - CITRUS PEEL FRAGRANCES
- EX03 - FRANKINCENSE FRAGRANCES
- EX04 - GARLIC FLAVORS
- EX05 - GARLIC FRAGRANCES
- EX06 - GINGER FLAVORS
- EX07 - GINGER FRAGRANCES
- EX08 - LAVENDER FLAVORS
- EX09 - LAVENDER FRAGRANCES
- EX10 - MINT FLAVORS
- EX11 - MINT FRAGRANCES
- EX12 - ORANGE PEEL FLAVORS
- EX13 - ORANGE PEEL FRAGRANCES
- EX14 - ROSE FRAGRANCES
- EX15 - ROSEMARY FLAVORS
- EX16 - ROSEMARY FRAGRANCES
- EX17 - TYMUS FRAGRANCES
- EX18 - STRAWBERRY FLAVORS

# Introduction

The microwave heating is generated by two mechanisms: ionic conduction and dipole rotation. The electromagnetic field applied by the microwave irradiation (MW) causes the migration of ions of the solvent and the resistance of this ion conduction results in friction that heats the solution. The electromagnetic field also produces a rearrangement of dipoles, which contributes to the heating of the solvent as well. This conversion of electromagnetic energy into thermal energy is peculiar and unique to this microwave heating process. In conventional processes the energy is transferred to the material by convection, conduction and radiation phenomena from the external material, while in microwave-assisted extraction (MAE) the microwave energy leads directly to molecular interactions with the electromagnetic field and consequent heating of the susceptible material<sup>[1]</sup>.

The MAE process occurs as the result of changes in the cell structure caused by electromagnetic waves. Its process acceleration and high extraction yield is the result of a synergistic combination of two transport phenomena: heat and mass gradients working in the same direction. In that sense MAE differs from the conventional extraction methods (solid-liquid or simply extraction) where the mass transfer occurs from the inside to the outside, although the heat transfer occurs from the outside to the inside of the substrate.



Microwave selective heating

Furthermore in MAE the dissipation of the heating takes place volumetrically inside the irradiated medium, while in conventional extraction the heat is transferred from the heating medium to the interior of the sample. There are six main steps that characterize the microwave extraction process conducted with a solvent: 1) penetration of the solvent into the solid matrix; 2) solubilization and/or break down of components; 3) transport of the solute out of the solid matrix by internal diffusion; 4) migration of the extracted solute from the external surface of the solid into the bulk solution by external diffusion; 5) movement of the extract with respect to the solid; 6) separation and discharge of the extract and solid<sup>[2]</sup>. Regarding the extraction process itself, three main subprocesses are known to take place: 1) desorption at a constant velocity of the substrate from the outer surface of the particle (equilibrium phase); 2) after an intermediary transition phase, a mass transfer by internal diffusion and convection takes place in the solid-liquid interface; 3) in the last step the solute must overcome the interactions that bind it to the matrix and diffuse into the extracting solvent (external diffusion). This last step is considered as the limiting step of the process, since the extraction rate is really low and due to its irreversibility<sup>[3]</sup>.

# Important parameters in Microwave-Assisted Extraction

## EFFECT OF EXTRACTION TIME

Extraction time represents another parameter to be taken into consideration during MAE. The period of heating is indeed an important factor in MAE, since a selective heating occurs and extraction times are very short (from a few minutes to a half-hour) in order to avoid thermal degradation or oxidation of the target compounds, that might be sensitive to overheating<sup>[4]</sup>. The overheating depends on the dielectric properties of the solvent and must be avoided especially in the case of thermolabile constituents. In that sense, when longer extraction time is required, multiple step extraction cycles are the best solution to enhance extraction yields, avoiding long heating<sup>[4b, 5]</sup>. In flavonoids extraction from *R. astragali* for instance, an increase in yield was found with exposure up to 25 minutes and then the extraction yield started to decrease<sup>[5]</sup>.

## EFFECT OF MICROWAVE POWER AND EXTRACTION TIME

Similar to the period of heating, microwave power and consequent temperature are factors that must be properly set in order to maximize the extraction efficiency, without causing the degradation of thermally sensitive compounds. In MAE temperature is controlled by incident microwave power that controls the amount of energy provided to the matrix, which is converted to heat energy in the dielectric material<sup>[2]</sup>. Microwave power is closely related to the extraction time required, the quantity and type of sample we want to obtain the extract from. It acts as a driving force because it causes localized heating in the plant matrix, which allows the solute to diffuse out from the plant cells and dissolve in the solvent. For that reason, an increase in the microwave power improves the extraction yield, resulting in a shorter extraction time<sup>[4b, 6]</sup>. On the other hand it must be taken into consideration that a high microwave power might cause a rapid breakage of the cell wall, ending up with the leaching out both of the target compounds and undesired impurities, in addition to causing the degradation of thermally sensitive compounds<sup>[7]</sup>.

## EFFECT OF THE SAMPLE CHARACTERISTICS

Another important parameter affecting the MAE process is the characteristics of the sample. It is well known the higher the contact surface area is, the higher the extraction efficiency results, since finer particles allow a deeper penetration of the microwave<sup>[8]</sup>. Prior grinding and homogenizing of the samples is therefore strongly recommended. In some cases soaking of the dried plant material in the extracting solvent before MAE has resulted in improved yield<sup>[7]</sup>. The moisture of the matrix acts as an extraction solvent, it evaporates with heat and causes internal ruptures of the cells, increasing the extraction yield<sup>[9]</sup>.

# Microwave Extraction Techniques

Due to the selective heating of microwaves accomplished by the microwave energy directly delivered to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy, the microwave extraction process can be defined in three specific extraction processes: Microwave Hydrodistillation (MWHD), Microwave Hydrodiffusion and Gravity (MHG) and Solvent-Free Microwave Extraction (SFME).

## MICROWAVE HYDRODISTILLATION (MWHD)

The microwave hydrodistillation system follows the classical hydrodistillation principle developed by Stashenko et al. in 2004<sup>[10]</sup>. The matrix is installed with water into a reactor that has already been placed inside the microwave oven. The refrigeration system and the part estimated to recover essences are situated outside the oven<sup>[2, 10]</sup>. This technique is widely applied for the extraction of volatile essential oils at atmospheric pressure from many aromatic plants and spices. During distillation, fragrance plants are exposed to boiling water, releasing their essential oils through evaporation. The recovery of essential oils is accomplished by the distillation apparatus placed outside the oven, based on the principle that the combined vapor pressure equals the ambient pressure at the boiling point, thus the ingredients of the essential oils can be evaporated at a temperature close to the water boiling point. After steam and essential oils have been condensed, both are collected in a so called Florentine flask where the essential oil floats on the top, while water goes to the bottom and can be easily separated<sup>[2]</sup>.


## SOLVENT-FREE MICROWAVE EXTRACTION (SFME)

Solvent-free microwave-assisted extraction (SFME) was developed and patented in 2004 by Chemat et al.<sup>[11] [12]</sup> Based on a relatively simple principle, this method involves placing fresh raw material in the microwave reactor, without any added solvent or water. The internal heating of the in situ water within the plant matrix distends the oil glands and sacs and leads to rupture of the glands and oleiferous receptacles. This process thus frees essential oil, which is evaporated by the in situ water of the plant material. A cooling system outside the microwave oven condenses the distillate continuously. Excess water is refluxed back to the extraction vessel in order to restore the water to the plant material, while essential oil is collected on the top of the microwave-clavenger module of the refrigerator system and can be easily collected out through the tap.

## MICROWAVE HYDRODIFFUSION AND GRAVITY (MHG)

Microwave hydrodiffusion and gravity (MHG) was patented by Chemat et al. in 2008<sup>[13]</sup>. This green extraction process is an original combination of microwave heating and gravity working at atmospheric pressure. MHG was conceived for the extraction of flavors from different plant materials at laboratory and industrial scales. In this relatively simple method, the plant material is directly placed in a microwave reactor without any added solvent or water. The internal heating of the in situ water within the plant material distends the plant cells and leads to rupture of the glands and cell receptacles. Heating under microwaves thus frees molecules of interest together with in situ water. This physical phenomenon, known as hydro-diffusion, allows the extract to diffuse outside the plant material and drop by earth gravity out of the microwave reactor through the perforated Pyrex disc. A heat exchanger outside the microwave oven cools the extract continuously. The extract (in situ water and metabolites) is collected and separated in receiving beakers. It is important to note this green method allows extracting of natural substances without distillation or solvent extraction which are the most energy- and solvent-consuming unit operations.

# How does Milestone conceive these extraction techniques?





**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX00 - BASIL FRAGRANCES

---

Basil Solvent-Free Microwave  
Extraction (SFME) and Microwave  
Hydrodistillation (MWHD)



### Introduction

*Ocimum Basilicum* L. (Lamiaceae), respectively, named basil, is an aromatic herb that has been used traditionally as a medicinal herb in the treatment of headaches, coughs, diarrhea, constipation, warts, worms and kidney malfunctions. It has a long history as culinary herb, thanks to its foliage adding a distinctive flavor to many foods. It is also a source of aroma compounds and essential oils containing biologically active constituents that possess insecticidal, nematocidal, fungistatic and antimicrobial properties<sup>[1]</sup>.

[1] O. Politeo, M. Jukic, M. Milos, *Food Chemistry* 2007, 101, 379-385.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from basil can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-Free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHD techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Fresh Basil, SFME						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1kW	2.1kW		
Small	500	500	•		0.4	0.07
Medium	1580	1580		•	1.3	0.08
Large	3720	1800		•	3	0.08

Dry Basil, MWHD						
Reactor	Weighted dry soaked material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1 kW	2.1kW		
Small	500	500	•		0.2	0.04
Medium	1580	1580		•	0.5	0.03
Large	3720	1800		•	1.5	0.04

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1



## Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows you to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

In this study, we propose state-of-the-art processes for the extraction of essential oils from basil through SFME

and MWHD. It is the unique modern concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

No.	Identified compound	Peak area (%)	RI <sup>a</sup> HP-20M	RI <sup>a</sup> HP-101
1	$\beta$ -Pinene	0.1	-	949
2	Limonene	0.1	1180	1005
3	1,8-Cineole	4.0	1185	1006
4	Camphor	0.5	1477	1109
5	Linalool	28.6	1518	1092
6	Bornyl acetate	0.5	1545	1252
7	Terpinen-4-ol	0.7	1563	1154
8	$\alpha$ -Bergamotene	2.2	1564	1407
9	Caryophyllene	0.3	-	1385
10	Alloaromadendrene	0.1	-	1450
11	Estragole	21.7	1632	1177
12	$\alpha$ -Terpineol	1.0	1653	1176
13	Germacrene D	0.3	1673	1444
14	$\alpha$ -Humulene	0.2	-	1417
15	Carvone	0.4	1685	1207
16	$\beta$ -Cubebene	0.5	1694	1059
17	$\beta$ -Burbonene	t	-	1354
18	$\beta$ -Elemene	0.3	-	1364
19	$\alpha$ -Cadinene	0.2	1716	1426
20	Calamenene	0.2	-	1483
21	$\alpha$ -Amorphene	1.0	1710	1479
22	$\beta$ -Farnesene	0.2	-	1452
23	$\Delta$ -Cadinene	0.1	1724	1486
24	$\alpha$ -Bisabolene	0.1	-	1506
25	(Z)-Methyl cinnamate	1.6	1900	1281

Table 2. Chemical composition of basil essential oil

No.	Identified compound	Peak area (%)	RI <sup>a</sup> HP-20M	RI <sup>a</sup> HP-101
26	Methyl eugenol	3.1	1959	1378
27	(E)-Methyl cinnamate	14.3	2019	1364
28	Spatulenol	0.8	2066	-
29	Eugenol	5.9	2105	1368
30	Carvacrol	t	2118	1814
31	■-Cadinol	7.1	2120	1614
32	Torreyol	0.2	2173	-
33	Chavicol	0.7	-	.. <sup>b</sup>
	Total	97.0		

–, not identified.

t, trace (<0.1%).

<sup>a</sup> Retention indices relative to C<sub>8</sub>–C<sub>22</sub> n-alkanes on polar HP-20M and apolar HP-101 column.

<sup>b</sup> Retention times is outside of retention times of homologous series of C<sub>8</sub>–C<sub>22</sub> n-alkanes (identified by MS).

*Table 2 (continued).*



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX01 - CITRUS PEEL FLAVOR

---

Citrus Peel Microwave  
Hydrodiffusion and Gravity (MHG)



### Introduction

Citrus is the most abundant crop in the world, with about 64 million tons of orange and 13 million tons of lemon products produced during 2004. The amount of residue obtained from citrus fruits account for 50% of the original amount of whole fruit. Produced in tones per day, citrus by-products represent a problem for management, pollution, and environmental issues, due to microbial spoilage. Citrus peel of fruit processing which provides a great potential for further commercial use. During the process of juice extraction oil sacs break and release volatile oils which are in pockets localized in the external part of the mesocarpe of fruit (flavedo). These oils are used in food and pharmaceutical industries, but can also provide flavoring ingredients to drinks, ice creams and other food products. In addition, substantial quantity of these oils are also used in the preparation of toilet soaps, perfumes, cosmetics and other home care products. Viro, Tomao, Ginies, Visinoni, and Chemat (2008) reported that the d-limonene, major component of the oil extracted from citrus peels, could be used as green solvent instead of hazardous petroleum

solvents for fats and oils determination. D-limonene is considered a very versatile chemical which can be used in a wide variety of applications <sup>[1][2]</sup>.

[1] N. Sahraoui, M. Abert Vian, M. El Maataoui, C. Boutekedjiret, F. Chemat, *Innovative Food Science and Emerging Technologies* 2011, 12, 163-170.

[2] V. Virost, V. Tomao, C. Ginies, F. Visinoni, F. Chemat, *Journal of Chromatography A* 2008, 1196, 147-152.

### Why Choose Microwave Flavor Set-up?

The patented and revolutionary Microwave Hydrodiffusion and Gravity (MHG) system paves the way to new flavoring products which were impossible to be obtained with the ancient extraction concepts. Percolation, solvent extraction, etc. were inefficient and not environmentally-friendly methods of flavor extraction. MHG improves the efficiency and the quality of flavoring products.

- New natural flavors
- Fast extraction
- No thermal degradation
- No solvent

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through MHG techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The MHG technique is suitable for both dry and fresh raw material, see the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Fresh Citrus Peel (MHG)							
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Total flavor extract [mL]	Volatile fraction [mL]	Total flavor extract yield [%]
			1kW	2.1kW			
Small	500	500	•		190	9.5	38
Medium	1580	1580		•	621	30	39.3
Large	3720	1800		•	1400	70.5	38.7

Dry Citrus Peel (MHG)							
Reactor	Weighted dry soaked material [g]	Power [W]	Chiller		Total flavor extract [mL]	Volatile fraction [mL]	Total flavor extract yield [%]
			1 kW	2.1kW			
Small	500	500	•		150	6.2	30
Medium	1580	1580		•	540	24	34.2
Large	3720	1800		•	1330	59	35.8

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g) for 20 min.

> 1800g: Power = 1800W for 40min

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

**Important Remarks**

Please take into account that the interaction between microwaves and raw material is based on the water content. The optimized method (power and time) depends on the type of citrus, mainly to the water content of citrus itself. Please use the reported parameter as general application note to start to optimize your own method. Be careful, using an excess power might cause burning of your sample. Please make sure to properly seal the glass reactor during the installation of the flavors set-up according to the manual, to avoid loss of vapor during extraction.

**Conclusion**

A newly and cleaner design process for extraction of flavors was developed in this study: MHG. This green process has been studied and tested using citrus peels. This new system was developed to indicate that the microwave extraction process of flavors offers important advantages over antiquated extraction techniques, namely, shorter extraction times, less energy consuming, lower costs as well as new flavoring products. Furthermore, the MHG system offers the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX02 - CITRUS PEEL FRAGRANCES

---

Citrus Peel Solvent-Free Microwave  
Extraction (SFME) and Microwave  
Hydrodistillation (MWHD)



### **Introduction**

Citrus essential oils are the most widely used essential oils in the world. They are obtained as by-products of the citrus processing. They are used as aroma flavor in many food products, including alcoholic and non-alcoholic beverages, candy and gelatins. In pharmaceutical industries they are employed as flavoring agents to mask unpleasant tastes of drugs. In perfumery and cosmetic, they are used in many preparations. The traditional way to extract essential oils is by cold pressing the citrus peels. The oil is present in oil sacs or oil glands located at different depths in the peel and the cuticles of the fruit. Peel and cuticle oils are removed mechanically by cold pressing and since cold pressing yields a watery emulsion, this emulsion is then centrifuged to separate out the essential oil. Distillation is also used in some countries as an economical way to recover the oils. During distillation, the citrus peel is exposed to boiling water or steam, releasing the essential oils through evaporation. Researchers in many universities are working on novel techniques that could lead to compact, safe, efficient,

energy saving, and sustainable extraction processes. Solvent-Free Microwave Extraction (SFME) as upcoming extraction techniques have been reported for the extraction of fragrances and flavors from citrus peel<sup>[3]</sup>.

[3] N. Bousbia, M. Abert Vian, M. Ferhat, B. Meklati, F. Chemat, Journal of Food Engineering 2009, 90, 409-413.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from citrus peel can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-Free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHD techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 2*).

Fresh Citrus Peel, SFME						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1kW	2.1kW		
Small	500	500	•		4.9	1.0
Medium	1580	1580		•	16	1.0
Large	3720	1800		•	48.4	1.3

Dry Citrus Peel, MWHD						
Reactor	Weighted dry soaked material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1 kW	2.1kW		
Small	500	500	•		2.8	0.56
Medium	1580	1580		•	9	0.57
Large	3720	1800		•	22.7	0.61

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

## Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

In this study, we propose state-of-the-art processes for the extraction of essential oils from citrus peel through SFME and MWHD. It is the unique modern concept of

the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

No.	Compounds <sup>a</sup>	R.I. <sup>b</sup>	R.I. <sup>c</sup>	SFME
	<b>Monoterpenes</b>			<b>92.76</b>
1	Pinene<Alpha->	926	1023	1.75
2	Pinene<Beta->	974	1109	15.35
3	Myrcene<Beta->	988	1165	1.33
4	Carene<Delta-3-> 1101	1101	1290	0.20
5	Limonene	1030	1206	65.25
6	Terpinene<Gamma->	1103	1285	8.08
	<b>Oxygenated Monoterpenes</b>			<b>92.76</b>
7	Linalool	1125	1538	0.18
8	Citron ellal	1167	1478	0.05
9	Terpin-4-ol	1191	1590	0.42
10	Terpineol<Alpha->	1203	1677	0.56
11	Nerol	1237	1781	0.49
12	Neral	1268	1670	0.68
13	Geraniol	1271	1828	0.60
14	Geranial	1284	1714	0.89
	<b>Sesquiterpenes</b>			<b>1.06</b>
15	Elemene<Beta->	1373	1583	-
16	Caryophellene<E->	1391	1594	0.18
17	Bergamotene<Alpha-Trans->	1437	1577	0.28
18	Humulene<Alpha->	1450	1657	0.04
19	Farnesene<(E)-Beta->	1453	1650	0.02
20	Germacrene D	1477	1696	-
21	Valencene	1488	1705	0.04
22	Bisabolene<(Z)-Alpha->	1498	1761	0.03
23	Bisabolene(Beta-)	1508	1718	0.44

Table 3. Chemical compositions of essential oils from citrus peel obtained by SFME



No.	Compounds <sup>a</sup>	R.I. <sup>b</sup>	R.I. <sup>c</sup>	SFME
	<b>Oxygenated Sesquiterpenes</b>			<b>0.03</b>
24	Elemol	1540	1381	-
25	Nerolidol<E->	1555	2026	-
25	Bisabolol<Alpha->	1684	2212	0.03
27	Nootkatone	1799	2250	-
	<b>Other oxygenated compounds</b>			<b>0.86</b>
28	Nonanal<N->	1126	1400	0.06
29	Citronellyl Acetate	1342	1645	0.04
30	Neryl Acetate	1351	1706	0.19
	Extraction time (min)			45
	Total oxygenated compounds (%)			4.78
	Total non oxygenated compounds (%)			93.82

<sup>a</sup> Essential oil compounds sorted by chemical families and percentages calculated by GC-FID on non-polar HP5MS™ capillary column.

<sup>b</sup> Retention indices calculated on non-polar HP5MS™ capillary column.

<sup>c</sup> Retention indices calculated on polar Carbowax™-PEG capillary column.

*Table 3 (continued).*



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX03 - FRANKINCENSE FRAGRANCES

### Hydrodistillation (MWHD)



#### Introduction

Gum resins from *Boswellia* species, also known as frankincense, have been used as a major ingredient in Ayurvedic and Chinese medicine to treat a variety of health-related conditions. Both frankincense chemical extracts and essential oil prepared from *Boswellia* species gum resins exhibit anti-neoplastic activity, and have been investigated as potential anti-cancer agents. The anti-cancer activity is mediated through multiple signaling pathways. In addition, frankincense essential oil overcomes multicellular resistant and invasive phenotypes of human breast cancer cells. Fast and green extraction of frankincense essential oil turns out to be extremely important. This essential oil is obtained from a resin from the bark of a shrub originally from the area surrounding the Red Sea, in Somalia and Arabia. To collect the resin, fine incisions are made in the bark, and drops of sap appear and dry in large, odorous yellow droplets<sup>[1]</sup>.

[1] X. Ni, M. Suhail, Q. Yang, A. Cao, K.-M. Fung, R. Postier, C. Woolley, G. Young, J. Zhang, H.-K. Lin, BMC Complementary and Alternative Medicine 2012, 12.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from frankincense can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) technique allows the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHD techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Frankincense, MWHD							
Reactor	Weighted dry material [g]	Weighted dry material + added water [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
				1 kW	2.1kW		
Small	600	1050	1050	•		27	4,5
Medium	1700	3318	1800		•	71,4	4,2
Large	4000	7816	1800		•	168	4,2

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

*Table 1*

## Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows you to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

In this study, we propose state-of-the-art processes for the extraction of essential oils from frankincense resin

through SFME and MWHD. It is the unique modern concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

RRI	Compounds	Boswellia rivae (%)			
1032	α-Pinene	5.3	1597	Bornyl acetate	0.4
1035	α-Thujene	1.3	1600	β-Elementene	0.3
1076	Camphene	0.1	1611	Terpinen-4-ol	1.4
1118	β-Pinene	0.6	1524	trans-Dihydrocarvone	0.2
1132	Sabinene	1.2	1639	trans-p-Mentha-2,8-diene-1-ol	3.9
1159	δ-3-Carene	9.6	1642	Thuj-3-en-10-al	0.2
1187	o-Cymene	2.5	1648	Myrtenal	0.5
1203	Limonene	14.8	1651	Sabinaketone	0.2
1213	1,8-Cineole	0.3	1657	Umbellulone	0.1
1266	(E)-β-Ocimene	0.4	1663	cis-Verbenol	0.5
1278	m-Cymene	0.4	1664	trans-Pinocarveol	2.2
1280	p-Cymene	2.9	1678	cis-p-Mentha-2,8-diene-1-ol	0.9
1424	o-Methylanisol	0.2	1683	trans-Verbenol	6.8
1430	α-Thujone	0.1	1700	p-Mentha-1,8-diene-4-ol	0.4
1439	γ-Campholene aldehyde	0.1	1706	α-Terpineol	1.4
1444	2,5-Dimethylstyrene	0.2	1709	α-Terpinyol acetate	1.0
1450	trans-Linalool oxide	0.1	1720	trans-Sabinol	0.3
1451	β-Thujone	0.7	1725	Verbenone	4.3
1458	cis-1,2-Limonene epoxide	4.6	1751	Carvone	1.6
1468	trans-1,2-Limonene epoxide	0.5	1804	Myrtenol	0.7
1474	trans-Sabinene hydrate	0.9	1811	trans-p-Mentha-1(7),8-diene-2-ol	0.3
1478	cis-Linalool oxide	0.1	1845	trans-Carveol	2.5
1498	(E)-β-Ocimene epoxide	0.2	1856	m-Cymen-8-ol	3.1
1499	α-Campholene aldehyde	0.9	1864	p-Cymen-8-ol	2.0
1536	Pinocamphone	0.3	1882	cis-Carveol	0.7
1553	Linalool	0.2	1896	cis-p-Mentha-1(7),8-diene-2-ol	0.2
1556	cis-Sabinene hydrate	0.9	1949	Piperitenone	0.6
1565	8,9-Limonene epoxide-I	0.6	2113	Cumin alcohol	0.1
1571	8,9-Limonene epoxide-II	0.6	2198	Thymol	0.1
1586	Pinocarvone	0.5	2239	Carvacrol	0.1
				Total	88.1

Table 2. The composition of *Boswellia Rivae* essential oil



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX04 - GARLIC FLAVORS

---

### Garlic Microwave Hydrodiffusion and Gravity (MHG)



#### **Introduction**

Garlic (*Allium Sativum*) belongs to a group of dietary supplements that lessen the incidence of cardiovascular and cerebrovascular diseases by reducing cholesterol concentration. The beneficial effect of garlic on health has been confirmed in several studies which showed that garlic has been evaluated for lowering blood pressure, cholesterol and glucose concentration, reducing blood lipids as well as for the prevention of arteriosclerosis and cancer. The biological activities of garlic including antibacterial, parasitocidal, antithrombotic, antioxidant and antidiabetic actions have been known for a long time. The unique flavor and health-promoting functions of garlic are generally attributed to its rich content of sulfur-containing compounds: alliin, g-glutamylcysteine and their derivatives.

### Why Choose Microwave Flavor Set-up?

The patented and revolutionary Microwave Hydrodiffusion and Gravity (MHG) system paves the way to new flavoring products which were impossible to be obtained with the ancient extraction concepts. Percolation, solvent extraction, etc. were inefficient and not environmentally-friendly methods of flavor extraction. MHG improves the efficiency and the quality of flavoring products.

- New natural flavors
- Fast extraction
- No thermal degradation
- No solvent

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the MHG techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The MHG technique is suitable for both dry and fresh raw material, see the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Fresh Garlic (MHG)						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Total flavor extract [mL]	Total flavor extract yield [%]
			1kW	2.1kW		
Small	1000	1000	•		185	18.5
Medium	3160	1800		•	600	19
Large	7445	1800		•	1400	18.8

Time, Power

≤ 1800 g: Power(W) = Weight(g) for 20 min.

> 1800g: Power = 1800W for 40min

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

### Important Remarks

Please take into account that the interaction between microwaves and raw material is based on the water content. The optimized method (power and time) depends on the type of garlic, mainly to the water content of garlic itself. Please use the reported parameter as a general application note to start to optimize your own method. Be careful, that using an excess power might cause burning of your sample. Please make sure to properly seal the glass reactor during the installation of the flavors set-up according to the manual, to avoid loss of vapor during extraction.

### Conclusion

A newly and cleaner design process for extraction of flavors was developed in this study: MHG. This green process has been studied and tested using garlic. This new system was developed to indicate that the microwave extraction process of flavors offers important advantages over antiquated extraction techniques, namely, shorter extraction times, less energy consuming, lower costs as well as new flavoring products. Furthermore, the MHG system offers the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX05 - GARLIC FRAGRANCES

---

Garlic Solvent-Free Microwave  
Extraction (SFME) and Microwave  
Hydrodistillation (MWHD)



### **Introduction**

Garlic (*Allium Sativum*) belongs to a group of dietary supplements that lessen the incidence of cardiovascular and cerebrovascular diseases by reducing cholesterol concentration. The beneficial effect of garlic on health has been confirmed in several studies which showed that garlic has been evaluated for lowering blood pressure, cholesterol and glucose concentration, reducing blood lipids as well as for the prevention of arteriosclerosis and cancer. The biological activities of garlic including antibacterial, parasitocidal, antithrombotic, antioxidant and antidiabetic actions have been known for a long time. The unique flavor and health-promoting functions of garlic are generally attributed to its rich content of sulfur-containing compounds: alliin, g-glutamylcysteine and their derivatives.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from garlic can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHd) and Solvent-Free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHd techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHd techniques are respectively suitable for fresh and dry raw materials. See the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 2*).

Fresh Garlic, SFME						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1kW	2.1kW		
Small	1000	1000	•		2	0.20
Medium	3160	1800		•	8.2	0.26
Large	7445	1800		•	22.3	0.30

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

*Table 2*



## Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows you to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

In this study, we propose state-of-the-art processes for the extraction of essential oils from garlic through SFME

and MWHD. It is the unique modern concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

Compounds	LRI <sub>HP1</sub>	LRI <sub>INNO</sub>	SFME (% ± SD)	Identification methods
Dimethyl sulfide	-	750	tr	SM, LRI
Dimethyl disulfide	731	1099	tr	SM, LRI, Std
Methyl ethyl disulfide	818	1119	0.5	SM, LRI
Diallyl sulfide	840	1150	0.5	SM, LRI, Std
Methyl allyl disulfide	894	1290	0.8 ± 0.1	SM, LRI, Std
(Z)-prop-1-enyl methyl disulfide	917	1251	tr	SM, LRI
(E)-prop-1-enyl methyl disulfide	919	1275	0.1	SM, LRI
Dimethyl trisulfide	941	1340	0.2	SM, LRI, Std
Diallyl disulfide	1056	1491	26.3 ± 0.9	SM, LRI, Std
Allyl (Z)-prop-1-enyl disulfide	1073	1390	2.9 ± 0.1	SM, LRI
Allyl (E)-prop-1-enyl disulfide	1082	1415	7.7 ± 0.1	SM, LRI
Allyl methyl trisulfide	1115	1601	7.9 ± 0.1	SM, LRI, Std
Methyl (E)-propenyl trisulfide	1138	-	0.1	Tentative
3-vinyl-(4H)-1,2-dithiin	1156	1735	1.1	SM, LRI
Unknown 1	1169	-	0.1	-
2-vinyl-(4H)-1,3-dithiin	1178	1830	3.1	SM, LRI
Diallyle trisulfide	1285	1825	39.7 ± 0.8	SM, LRI, Std
Propyl propenyl trisulfide <sup>d</sup>	1290	1781	0.3	-
Allyl propenyl trisulfide	1300	1798	2.0	SM, LRI
3,5-diethyl-1,2,4-trithiolane	1321	1788	0.2	SM, LRI
Unknown 2	1348	2011	0.4	-
Diallyl tetrasulfide	1508	-	1.5	SM, LRI, Std
2,4-dimethyl-5,6-dithia-2,7-nonadienal	1718	> 2400	0.9 ± 0.1	SM, LRI
2,4-dimethyl-5,6-dithia-2,7-nonadienal	1730	> 2400	0.4 ± 0.1	SM, LRI
Unknown 3	1775	-	0.1	-

<sup>a</sup> Compounds are listed in order of their elution time from a HP-1 column. Compositional values less than 0.1% are denoted as traces (tr). Presence of a compound is indicated by its GC-FID percentage with S.D., absence is indicated by "-".

<sup>b</sup> RI = retention indices are determined on HP-1 and INNOWAX column using the homologous series of n-alkanes (C6-C24).

<sup>c</sup> S.D. = standard deviation.

<sup>d</sup> Correct isomer not identified

tentative: tentatively identified by MS and RI without standard compound co-injection

Unknown 1: 186 (M+, 7.1); 162(23.2); 100 (15); 97 (28.3); 73(50.1); 60 (16.0); 56(15.6); 59 (28.2); 57 (31.8); 45 (45.5); 41(100).

Unknown 2: 178 (M+, 3.4); 172(16.6); 170 (100); 128(42.3); 106(25.3); 73(12.3); 64(62.3); 59(17.3); 45(26.8); 42 (11.7); 41(31).

Unknown 3: 179(M+, 5.6); 147(89.0); 106(10.5); 105(33.4); 75(14.6); 73(98.7); 64 (14.8); 57(10.7); 47 (15.4); 45(29.7); 41(100).

Table 3. Essential oil from garlic (*Allium sativum* oil) obtained by SFME



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX06 - GINGER FLAVORS

### Ginger Microwave Hydrodiffusion and Gravity (MHG)



#### Introduction

Due to its composition in valuable natural compounds, ginger represents a suitable matrix for extraction of essential oil as well as many other bioactive compounds. It contains products of interest such as essential oils (1–4%), phenolics (gingerols and 6-shogaol, 1–2%), and total carbohydrates (60–75%). Ginger, and more specifically rhizomes, are variously used as food product or in traditional medicine. In the food industry, rhizomes are mainly used for spices or condiments (fresh or dried), candy or as juice after cold mechanical pressing. Due to the fact that mechanical pressing does not alter the chemical composition of the pressed product, this process provides huge amounts of press cake still containing high amounts of bioactive compounds, which is currently considered as waste. The novelty of the extraction of these products through MHG and SFME relies on the extraction of compounds achieved without addition of solvent or water. The only water used in the process is the matrix water naturally present in the plant cells <sup>[1]</sup>.

### Why Choose Microwave Flavor Set-up?

The patented and revolutionary Microwave Hydrodiffusion and Gravity (MHG) system paves the way to new flavoring products which were impossible to be obtained with the ancient extraction concepts. Percolation, solvent extraction, etc. were inefficient and not environmentally-friendly methods of flavor extraction. MHG improves the efficiency and the quality of flavoring products.

- New natural flavors
- Fast extraction
- No thermal degradation
- No solvent

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the MHG techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The MHG technique is suitable for both dry and fresh raw material, see the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Fresh Ginger (MHG)						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Total flavor extract [mL]	Total flavor extract yield [%]
			1kW	2.1kW		
Small	1000	1000	•		388	38.8
Medium	3160	1800		•	1230	39
Large	7445	1800		•	2910	39.1

Time, Power

≤ 1800 g: Power(W) = Weight(g) for 20 min.

> 1800g: Power = 1800W for 40min

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

### Important Remarks

Please take into account that the interaction between microwaves and raw material is based on the water content. The optimized method (power and time) depends on the type of ginger, mainly to the water content of ginger itself. Please use the reported parameter as a general application note to start to optimize your own method. Be careful, that using an excess power might cause burning of your sample. Please make sure to properly seal the glass reactor during the installation of the flavors set-up according to the manual, to avoid loss of vapor during extraction.

### Conclusion

A newly and cleaner design process for extraction of flavors was developed in this study: MHG. This green process has been studied and tested using ginger. This new system was developed to indicate that the microwave extraction process of flavors offers important advantages over antiquated extraction techniques, namely, shorter extraction times, less energy consuming, lower costs as well as new flavoring products. Furthermore, the MHG system offers the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX07 - GINGER FRAGRANCES

---

Ginger Solvent-Free Microwave  
Extraction (SFME) and Microwave  
Hydrodistillation (MWHD)



### Introduction

Due to its composition in valuable natural compounds, ginger represents a suitable matrix for extraction of essential oil as well as many other bioactive compounds. It contains products of interest such as essential oils (1–4%), phenolics (gingerols and 6-shogaol, 1–2%), and total carbohydrates (60–75%). Ginger, and more specifically rhizomes, are variously used as food product or in traditional medicine. In the food industry, rhizomes are mainly used for spices or condiments (fresh or dried), candy or as juice after cold mechanical pressing. Due to the fact that mechanical pressing does not alter the chemical composition of the pressed product, this process provides huge amounts of press cake still containing high amounts of bioactive compounds, which is currently considered as waste. The novelty of the extraction of these products through MHG and SFME relies on the extraction of compounds achieved without addition of solvent or water. The only water used in the process is the matrix water naturally present in the plant cells <sup>[1]</sup>.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from ginger can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHHD) and Solvent-Free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHHD techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHHD techniques are respectively suitable for fresh and dry raw materials. See the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 2*).

Fresh Ginger, SFME						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1kW	2.1kW		
Small	1000	1000	•		0.7	0.08
Medium	3160	1800		•	2.2	0.07
Large	7445	1800		•	5.2	0.07

Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

## Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows you to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

In this study, we propose state-of-the-art processes for the extraction of essential oils from ginger through SFME

and MWHD. It is the unique modern concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

	GR	GP	0.6 W g <sup>-1</sup>	0.8 W g <sup>-1</sup>	1.0 W g <sup>-1</sup>	1.2 W g <sup>-1</sup>	1.4 W g <sup>-1</sup>	1.6 W g <sup>-1</sup>	1.8 W g <sup>-1</sup>
Essential oil Yield (g per 100 g fresh plant material). Major compounds (%)	0.2	0.3	0.2	0.2	0.2	0.2	0.2	0.2	0.2
β-Pinene	1.2	1.0	2.3	2.6	2.4	2.6	2.3	2.4	2.2
Camphene	4.3	3.8	9.1	10.3	9.2	10.0	9.1	9.4	9.1
Sabinene	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Sulcatone	0.0	0.8	1.2	2.8	3.3	3.2	3.0	3.2	2.9
Myrcene	0.6	0.6	0.0	1.4	1.4	1.4	1.3	1.3	1.1
■-Phellandrene	0.2	0.1	0.2	0.3	0.3	0.3	0.3	0.3	0.1
Limonene	0.9	0.9	1.7	1.9	1.9	1.9	1.7	1.8	1.7
β-Phellandrene	4.6	4.2	8.7	10.4	10.3	10.2	9.7	10.0	8.6
Terpinolene	0.1	0.1	0.1	0.1	0.3	0.3	0.3	0.3	0.2
Linalol	0.2	0.2	0.3	0.4	0.4	0.4	0.4	0.4	0.4
Borneol	0.5	0.6	0.8	0.9	1.0	0.9	1.0	1.0	1.1
■-Terpineol	0.2	0.3	0.4	0.5	0.5	0.5	0.5	0.5	0.6
Citronellol	0.1	0.3	0.2	0.5	0.4	0.3	0.4	0.4	0.8
Neral	1.7	0.5	0.4	1.3	1.5	1.7	1.5	1.5	1.3
Geraniol	0.1	0.2	0.1	0.3	0.3	0.2	0.2	0.2	0.6
Geranial	3.3	1.0	0.6	1.9	2.2	2.6	2.3	2.5	2.3
Geranyl acetate	0.3	0.1	0.4	0.2	0.2	0.2	0.2	0.2	0.2
■-Curcumene	3.5	13.9	17.0	7.6	7.2	6.6	7.0	6.8	9.9
Germacrene D	1.6	1.3	0.1	1.3	1.4	1.4	1.4	1.4	0.7
Zingiberene	35.7	25.2	18.4	23.2	24.0	24.0	25.1	24.3	18.4
■-Farnesene	6.5	6.5	6.3	5.4	5.5	5.5	5.7	5.5	5.7
β-Bisabolene	5.7	6.8	0.0	4.8	4.7	4.6	4.8	4.7	5.4
β-Sesquiphellandrene	12.1	13.9	12.3	9.9	9.9	9.7	10.2	9.8	10.4
Antioxidants Total content (g per 100 g plant material DW). Major compounds (g per 100 g plant material DW)	1.17	0.90	0.57	1.24	1.06	1.18	1.22	1.37	1.18
6-Gingerol	0.77	0.58	0.31	0.81	0.65	0.79	0.81	0.92	0.79
8-Gingerol	0.15	0.11	0.07	0.14	0.11	0.14	0.14	0.17	0.14
10-Gingerol	0.23	0.19	0.11	0.18	0.19	0.19	0.19	0.21	0.19
6-Shogaol	0.02	0.02	0.08	0.11	0.10	0.08	0.09	0.08	0.08

Table 3. Volatile compounds and antioxidants extracted from ginger plant material by SFME



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX08 - LAVENDER FLAVORS

---

### Lavender Microwave Hydrodiffusion and Gravity (MHG)



#### Introduction

Lavender is one of the most useful medicinal plants. Commercially, lavender provides several important essential oils to the fragrance industry, including soaps, colognes, perfumes, skin lotions and other cosmetics. In food manufacturing, lavender essential oil is employed in flavoring beverages, ice cream, candy, baked goods and chewing gum<sup>[1]</sup>. The essential oils of *Lavandula* species are obtained by steam distillation of the fresh flowering spikes. Oil quality is assessed by oil chemical composition and by the organoleptic opinion. In addition, a large range of medical uses for this plant have also been reported. These include antispasmodic, sedative, antihypertensive, antiseptic, healing and anti-inflammatory properties, all of which render it highly appreciated in phytotherapy and aromatherapy.

[1] F. Chemat, M. Lucchesi, J. Smadja, L. Favretto, G. Colnaghi, F. Visinoni, *Analytica Chimica Acta* 2006, 555, 157-160.

## Why Choose Microwave Flavor Set-up?

The patented and revolutionary Microwave Hydrodiffusion and Gravity (MHG) system paves the way to new flavoring products which were impossible to be obtained with the ancient extraction concepts. Percolation, solvent extraction, etc. were inefficient and not environmentally-friendly methods of flavor extraction. MHG improves the efficiency and the quality of flavoring products.

- New natural flavors
- Fast extraction
- No thermal degradation
- No solvent

## Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the MHG techniques (see the "Microwave Extraction Techniques" section for theory and principle).

## Results and experimental procedure

The MHG technique is suitable for both dry and fresh raw material, see the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Fresh Lavender (MHG)							
Reactor	Weighted fresh raw material [g]*	Power [W]	Chiller		Total flavor extract [mL]	Volatile fraction [mL]	Total flavor extract yield [%]
			1kW	2.1kW			
Small	1000	1000	•		360	6	36
Medium	3160	1800		•	950	20	30.1
Large	7445	1800		•	2500	39	33.6

Dry Lavender (MHG)							
Reactor	Weighted fresh raw material [g]*	Power [W]	Chiller		Total flavor extract [mL]	Volatile fraction [mL]	Total flavor extract yield [%]
			1 kW	2.1kW			
Small	1000	1000	•		285	4.8	28.5
Medium	3160	1800		•	882	15.2	27.9
Large	7445	1800		•	2144	35.7	28.8

\*Time, Power

≤ 1800 g: Power(W) = Weight(g) for 20 min.

> 1800g: Power = 1800W for 40min

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

## Important remarks

Please take into account that the interaction between microwaves and raw material is based on the water content. The optimized method (power and time) depends on the type of lavender, mainly to the water content of lavender itself. Please use the reported parameter as a general application note to start to optimize your own method. Be careful, that using an excess power might cause burning of your sample. Please make sure to properly seal the glass reactor during the installation of the flavors set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

A newly and cleaner design process for extraction of flavors was developed in this study: MHG. This green process has been studied and tested using lavender flowers. This new system was developed to indicate that the microwave extraction process of flavors offers important advantages over antiquated extraction techniques, namely, shorter extraction times, less energy consuming, lower costs as well as new flavoring products. Furthermore, the MHG system offers the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.





**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX09 - LAVENDER FRAGRANCES

---

### Lavender Solvent-Free Microwave Extraction (SFME) and Microwave Hydrodistillation (MWHD)



#### **Introduction**

Lavender is one of the most useful medicinal plants. Commercially, lavender provides several important essential oils to the fragrance industry, including soaps, colognes, perfumes, skin lotions and other cosmetics. In food manufacturing, lavender essential oil is employed in flavoring beverages, ice cream, candy, baked goods and chewing gum<sup>[1]</sup>. The essential oils of *Lavandula* species are obtained by steam distillation of the fresh flowering spikes. Oil quality is assessed by oil chemical composition and by the organoleptic opinion. In addition, a large range of medical uses for this plant have also been reported. These include antispasmodic, sedative, antihypertensive, antiseptic, healing and anti-inflammatory properties, all of which render it highly appreciated in phytotherapy and aromatherapy.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from lavender can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHHD) and Solvent-Free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHHD techniques (see the “Microwave Extraction Techniques” section for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHHD techniques are respectively suitable for fresh and dry raw materials. See the “Quick start guide” for a list of easy and sequential setting-up operations (*Table 2*).

Fresh Lavender, SFME						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1kW	2.1kW		
Small	1000	1000	•		2.5	0.25
Medium	3160	1800		•	6.7	0.21
Large	7445	1800		•	15	0.22

Dry Lavender, MWHHD						
Reactor	Weighted dry soaked material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1 kW	2.1kW		
Small	1000	1000	•		1.7	0.17
Medium	3160	1800		•	6.6	0.21
Large	7445	1800		•	14.2	0.19

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

## Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows you to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the flavors set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

In this study, we propose state-of-the-art processes for the extraction of essential oils from lavender through

SFME and MWHD. It is the unique modern concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

No.	Compounds	R.R.I. <sup>a</sup>	MWHD <sup>b</sup> (%)
	Monoterpenes		3.54
1	α-Thujene	907	0.08
2	α-Pinene	914	0.51
3	Camphene	933	0.32
4	Sabinene	963	0.14
5	β-Pinene	968	0.59
6	β-Myrcene	985	0.50
7	3-Carene	1008	0.22
8	Limonene	1024	tr.
9	(Z)-β-Ocimene	1031	0.33
10	(E)-β-Ocimene	1040	0.37
11	γ-Terpinene	1050	0.09
12	Terpinolene	1078	0.37
	Oxygenated monoterpenes		78.29
13	1,8-Cineole	1027	7.23
14	Sabinene hydrate-cis	1058	0.66
15	Linalool	1099	47.82
16	Camphor	1137	11.82
17	Borneol	1161	4.15
18	Terpin-4-ol	1174	5.94
19	p-Cymen-8-ol	1179	tr.
20	α-Terpineol	1186	0.68
	Sesquiterpenes		2.77
21	α-Bergamotene cis	1400	0.10
22	β-Caryophyllene	1412	1.28

Table 3. Yields, extraction times, grouped compounds and chemical compositions of essential oils obtained by MWHD from lavender flowers.

No.	Compounds	R.R.I. <sup>a</sup>	MWHD <sup>b</sup> (%)
23	■-Santalene	1414	0.15
24	(E)-β-Farnesene	1453	0.63
25	Sesquiterpene 1	1474	0.61
	Oxygenated sesquiterpenes		0.29
26	Caryophyllene oxide	1573	0.11
27	■-Bisabolol	1677	0.18
	Other oxygenated compounds		15.01
28	Octan-3-one	977	0.78
29	Octan-3-ol	990	0.26
30	Dihydromyrcenol	1063	0.34
31	n.i.	1141	0.37
32	n.i.	1188	2.00
33	n.i.	1232	0.43
34	Linalool acetate	1254	10.74
35	Geranyl acetate	1377	0.08
	Total extraction time (min)		45

tr., trace; n.i., non-identified.

<sup>a</sup> R.R.I., relative retention indices relative to C<sub>8</sub>–C<sub>22</sub> n-alkanes on SBP5™ capillary column.

<sup>b</sup> MWHD, microwave hydrodistillation.

Table 3 (continued).



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX10 - MINT FLAVORS

---

### Mint Microwave Hydrodiffusion and Gravity (MHG)



#### **Introduction**

The mint species have a great importance, both medicinal and commercial. Indeed, leaves, flowers and stems of *Mentha* spp. are frequently used in herbal teas or as additives in commercial spice mixtures for many foods to offer aroma and flavor. In addition, *Mentha* spp. has been used as a folk remedy for treatment of nausea, bronchitis, flatulence, anorexia, ulcerative colitis, and liver complaints due to its anti-inflammatory, carminative, antiemetic, diaphoretic, antispasmodic, analgesic, stimulant, emmenagogue, and anticephalalgic activities. Commercially, the most important mint species are peppermint (*M. x piperita*), spearmint (*M. spicata*), and corn mint (*M. canadensis*). From these species, corn mint is cultivated only because of oil production. Peppermint oil is one of the most popular and widely used essential oils, mostly because of its main components menthol and menthone. Corn mint is the richest source of natural menthol (Sharma and Tyagi 1991; Shasany et al. 2000). Carvone-scented mint plants, such as spearmint, are rich in carvone and are widely used as spices and cultivated in several countries. Peppermint oil is used for

flavoring pharmaceuticals and oral preparations, such as toothpastes, dental creams, and mouth washes. It is also used as a flavoring agent in cough drops, chewing gums, confectionery and alcoholic liqueurs. It is used in medicines for internal use. Its pleasant taste makes it an excellent gastric stimulant [1].

[1] H. Hajjaoui, N. Trabelsi, E. Noumi, M. Snoussi, H. Fallah, R. Ksouri, A. Bakhrouf, World Journal of Microbiology and Biotechnology 2009, 25, 2227-2238.

### Why Choose Microwave Flavor Set-up?

The patented and revolutionary Microwave Hydrodiffusion and Gravity (MHG) system paves the way to new flavoring products which were impossible to be obtained with the ancient extraction concepts. Percolation, solvent extraction, etc. were inefficient and not environmentally-friendly methods of flavor extraction. MHG improves the efficiency and the quality of flavoring products.

- New natural flavors
- Fast extraction
- No thermal degradation
- No solvent

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the MHG techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The MHG technique is suitable for both dry and fresh raw material, see the "Quick start guide" for a list of easy and sequential setting-up operations (Table 1).

Fresh Mint (MHG)						
Reactor	Weighted fresh raw material [g]*	Power [W]	Chiller		Total flavor extract [mL]	Total flavor extract yield [%]
			1kW	2.1kW		
Small	500	500	•		125	25
Medium	1580	1580		•	382	24.2
Large	3720	1800		•	960	25.8

Dry Mint (MHG)						
Reactor	Weighted fresh dry soaked material [g]*	Power [W]	Chiller		Total flavor extract [mL]	Total flavor extract yield [%]
			1 kW	2.1kW		
Small	500	500	•		110	22
Medium	1580	1580		•	302	19.1
Large	3720	1800		•	810	21.8

Time, Power

≤ 1800 g: Power(W) = Weight(g) for 20 min.

> 1800g: Power = 1800W for 40min

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

### **Important Remarks**

Please take into account that the interaction between microwaves and raw material is based on the water content. The optimized method (power and time) depends on the type of mint, mainly to the water content of mint itself. Please use the reported parameter as a general application note to start to optimize your own method. Be careful, that using an excess power might cause burning of your sample. Please make sure to properly seal the glass reactor during the installation of the flavors set-up according to the manual, to avoid loss of vapor during extraction.

### **Conclusion**

A newly and cleaner design process for extraction of flavors was developed in this study: MHG. This green process has been studied and tested using mint. This new system was developed to indicate that the microwave extraction process of flavors offers important advantages over antiquated extraction techniques, namely, shorter extraction times, less energy consuming, lower costs as well as new flavoring products. Furthermore, the MHG system offers the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX11 - MINT FRAGRANCES

---

Mint Solvent-Free Microwave  
Extraction (SFME) and Microwave  
Hydrodistillation (MWHD)



### Introduction

The mint species have a great importance, both medicinal and commercial. Indeed, leaves, flowers and stems of *Mentha* spp. are frequently used in herbal teas or as additives in commercial spice mixtures for many foods to offer aroma and flavor. In addition, *Mentha* spp. has been used as a folk remedy for treatment of nausea, bronchitis, flatulence, anorexia, ulcerative colitis, and liver complaints due to its anti-inflammatory, carminative, antiemetic, diaphoretic, antispasmodic, analgesic, stimulant, emmenagogue, and anticatharrhal activities. Commercially, the most important mint species are peppermint (*M. x piperita*), spearmint (*M. spicata*), and corn mint (*M. canadensis*). From these species, corn mint is cultivated only because of oil production. Peppermint oil is one of the most popular and widely used essential oils, mostly because of its main components menthol and menthone. Corn mint is the richest source of natural menthol (Sharma and Tyagi 1991; Shasany et al. 2000). Carvone-scented mint plants, such as spearmint, are rich in carvone and are widely used as spices and cultivated in several countries. Peppermint oil is used for



cultivated in several countries. Peppermint oil is used for flavoring pharmaceuticals and oral preparations, such as toothpastes, dental creams, and mouth washes. It is also used as a flavoring agent in cough drops, chewing gums, confectionery and alcoholic liqueurs. It is used in medicines for internal use. Its pleasant taste makes it an excellent gastric stimulant [1].

[1] H. Hajlaoui, N. Trabelsi, E. Noumi, M. Snoussi, H. Fallah, R. Ksouri, A. Bakhrouf, World Journal of Microbiology and Biotechnology 2009, 25, 2227-2238.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from mint can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-Free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHD techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the "Quick start guide" for a list of easy and sequential setting-up operations (Table 2).

Fresh Mint, SFME						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1kW	2.1kW		
Small	500	500	•		4.8	0.98
Medium	1580	1580		•	16	1.01
Large	3720	1800		•	37.2	1.0

Dry Mint, MWHD						
Reactor	Weighted dry soaked material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1 kW	2.1kW		
Small	500	500	•		4.3	0.86
Medium	1580	1580		•	13.7	0.87
Large	3720	1800		•	30.9	0.83

Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

## Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows you to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

In this study, we propose state-of-the-art processes for the extraction of essential oils from mint through SFME

and MWHD. It is the unique modern concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large) complying with a high range of extraction-scale needs.

	Fresh leaves SFME (g Kg <sup>-1</sup> )	Dried leaves MWHD (g Kg <sup>-1</sup> )
α-Pinene	4.0	7.4
Sabinene	4.6	6.3
β-Pinene	7.2	10.8
2-Thujene	15.9	14.7
3-Octanol	-	2.6
Limonene	67.5	8.67
1,8-Cineole	24.8	3.45
(E)-β-Ocimene	6.3	5.2
(Z)-β-Ocimene	3.4	2.0
γ-Terpinene	3.6	1.1
3-Carene	8.1	11.1
p-Menth-1-en-8-ol	-	1.4
(-)-4-Terpineol	7.1	2.6
t-Dihydrocarvone	18.9	17.3
c-Carveol	7.1	3.3
t-Carveol	-	2.3
D-Carvone	602.7	601.7
Piperitone	4.9	8.2
t-Carvone oxide	2.8	4.0
c-Carvone oxide	-	4.2
Dihydroedulan II	-	1.0
Dihydroedulan I	-	1.5
Isolimonene	5.8	-
Dihydrocarvyl acetate	-	3.2
t-Carveyl acetate	-	1.8
β-Bourbonene	29.7	24.7
β-Elemene	4.2	5.5
Isocaryophyllene	-	1.1

Table 3. Main components of *M. spicata* L. var. *rubra* EO yield (g kg<sup>-1</sup>)

	Fresh leaves SFME (g Kg <sup>-1</sup> )	Dried leaves MWHD (g Kg <sup>-1</sup> )
β-Caryophyllene	44.8	42.3
■-Caryophyllene	3.9	3.8
(+)-Epi-bicyclosesquiphellandrene	9.8	9.0
Germacrene D	17.3	16.1
Bicyclogermacrene	4.4	3.5
■-Muurolene	-	1.3
Calamenene	9.3	5.6
Caryophyllene oxide	16.8	3.0
1,4-Cadinadiene	5.9	2.2
τ -Muurolol	4.1	1.2
Total	945.0	954.0

*Table 3 (continued).*



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX12 - ORANGE PEEL FLAVORS

---

Orange Peel Microwave  
Hydrodiffusion and Gravity (MHG)



### Introduction

Orange essential oil is used to confer the aroma and flavor of orange to a wide variety of products, such as carbonated drinks, ice-creams, cakes, air-fresheners, perfumes and more. They are also being used in the design of new products, to which they add aroma and flavor 2, 3 and 4. Another application takes advantage of the germicide properties of some of their components. In that sense, a small amount of d-limonene was very effective in the germicide treatment of waste waters. Carotenoid pigments, found in orange extracts, are important for health, not only because of their nutritional value as precursors of vitamin A, but also because of their antioxidant potential. They also seem to have certain anticarcinogenic properties. These compounds are used in food coloring <sup>[1]</sup>.

[1] B. Mira, M. Blasco, A. Berna, S. Subirats, Journal of Supercritical Fluids 1999,

## Why Choose Microwave Flavor Set-up?

The patented and revolutionary Microwave Hydrodiffusion and Gravity (MHG) system paves the way to new flavoring products which were impossible to be obtained with the ancient extraction concepts. Percolation, solvent extraction, etc. were inefficient and not environmentally-friendly methods of flavor extraction. MHG improves the efficiency and the quality of flavoring products.

- New natural flavors
- Fast extraction
- No thermal degradation
- No solvent

## Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the MHG techniques (see the "Microwave Extraction Techniques" section for theory and principle).

## Results and Experimental Procedure

The MHG technique is suitable for both dry and fresh raw material, see the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Fresh Orange peel (MHG)							
Reactor	Weighted fresh raw material [g] <sup>*</sup>	Power [W]	Chiller		Total flavor extract [mL]	Volatile fraction [mL]	Total flavor extract yield [%]
			1kW	2.1kW			
Small	500	500	•		185	9.5	37
Medium	1580	1580		•	599	30	38
Large	3720	1800		•	1432	70.5	38.5

Dry Orange peel (MHG)							
Reactor	Weighted dry soaked material [g] <sup>*</sup>	Power [W]	Chiller		Total flavor extract [mL]	Volatile fraction [mL]	Total flavor extract yield [%]
			1 kW	2.1kW			
Small	500	500	•		150	6.2	30
Medium	1580	1580		•	540	24	34.2
Large	3720	1800		•	1330	59	35.8

<sup>\*</sup>Time, Power

≤ 1800 g: Power(W) = Weight(g) for 20 min.

> 1800g: Power = 1800W for 40min

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

## Important Remarks

Please take into account that the interaction between microwaves and raw material is based on the water content. The optimized method (power and time) depends on the type of orange fruit, mainly to the water content of orange fruit itself. Please use the reported parameter as general application note to start to optimize your own method. Be careful, that using an excess power might cause burning of your sample. Please make sure to properly seal the glass reactor during the installation of the flavors set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

A newly and cleaner design process for extraction of flavors was developed in this study: MHG. This green process has been studied and tested using orange peel. This new system was developed to indicate that the microwave extraction process of flavors offers important advantages over antiquated extraction techniques, namely, shorter extraction times, less energy consuming, lower costs as well as new flavoring products. Furthermore, the MHG system offers the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX13 - ORANGE PEEL FRAGRANCES

---

Orange Peel Solvent-Free  
Microwave Extraction (SFME) and  
Microwave Hydrodistillation (MWHD)



### Introduction

Orange essential oil is used to confer the aroma and flavor of orange to a wide variety of products, such as carbonated drinks, ice-creams, cakes, air-fresheners, perfumes and more. They are also being used in the design of new products, to which they add aroma and flavor 2, 3 and 4. Another application takes advantage of the germicide properties of some of their components. In that sense, a small amount of d-limonene was very effective in the germicide treatment of waste waters. Carotenoid pigments, found in orange extracts, are important for health, not only because of their nutritional value as precursors of vitamin A, but also because of their antioxidant potential. They also seem to have certain anticarcinogenic properties. These compounds are used in food coloring <sup>[1]</sup>.

[1] B. Mira, M. Blasco, A. Berna, S. Subirats, Journal of Supercritical Fluids 1999, 14, 95-104.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from orange peel can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHHD) and Solvent-Free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHHD techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHHD techniques are respectively suitable for fresh and dry raw materials. See the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 2*).

Fresh Orange peel, SFME						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1kW	2.1kW		
Small	500	500	•		4.9	1.0
Medium	1580	1580		•	16	1.0
Large	3720	1800		•	69	1.9

Dry Orange peel, MWHHD						
Reactor	Weighted dry soaked material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1 kW	2.1kW		
Small	500	500	•		2.8	0.56
Medium	1580	1580		•	9	0.57
Large	3720	1800		•	23.1	0.62

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

### Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows you to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

### Conclusion

In this study, we propose state-of-the-art processes for the extraction of essential oils from orange peel

through SFME and MWHD. It is the unique modern concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

N°	Compound	I <sup>a</sup>	SFME (%)
Monoterpenes			
1	$\alpha$ -Pinene	928	0.60
2	Sabinene	968	0.23
3	$\beta$ -Myrcene	988	1.81
4	3-Carene	1001	tr
5	$\delta$ -3-Carene	1007	0.05
6	Limonene	1034	96.20
7	$\alpha$ -Terpinolene	1087	0.01
Oxygenated monoterpenes			
8	Linalool	1093	0.17
9	Trans-Limonene oxide	1135	0.02
10	Citronellal	1150	0.02
11	$\beta$ -Citronellol	1227	0.06
Sesquiterpene			
12	$\alpha$ -Copaene	1372	0.01
13	$\beta$ -Cubebene	1382	tr
14	$\beta$ -elemene	1386	0.01
15	Caryophyllene (E)	1417	0.01
16	$\alpha$ -Humulene	1452	0.05
17	Germacrene-D	1479	tr
18	Valencene	1490	0.01
19	Germacrene-A	1503	0.01
20	$\delta$ -Cadinene	1520	0.01

Table 3. Chemical composition of orange EOs obtained by SFME



N°	Compound	I <sup>a</sup>	SFME (%)
Oxygenated sesquiterpenes			
21	Caryophyllene oxide	1589	tr
22	Cis, trans -Farnesol	1694	-
23	$\alpha$ -Sinensal	1754	tr
Other oxygenated compounds			
24	Decanal	1203	0.16
25	n-Dodecanal	1404	0.02
Extraction time (min.)			45

tr. tracesb0.01%.

<sup>a</sup> I, Retention indices relative to C5–C28 n-alkanes calculated on non-polar HP5MS capillary column. Percentages calculated by GC–FID on non-polar HP5MS capillary column. Essential oil compounds sorted by chemical families.

*Table 3 (continued).*



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX14 - ROSE FRAGRANCES

---

Rosa Damascena Mill. Solvent-Free  
Microwave Extraction (SFME) and  
Microwave Hydrodistillation (MWHD)



### **Introduction**

Roses have been used since the earliest times in rituals, cosmetics, perfumes, medicines and aromatherapy. A great variety of garden roses also exist, which are bred less for fragrance and more for color and shape. Even with the high price of roses and the advent of organic synthesis, rose oils are still the most widely used essential oils in perfumery. In fact, because of the labor-intensive production process and the low content of oil in rose blooms, rose oil is very expensive and is often called 'liquid gold'. For the production of essential rose oil as well as rose extracts, the two rose species most often used are Rosa Damascena (pink damask rose) and Rosa Centifolia (light pink cottage rose). Whereas the former is predominantly used for rose oil production, the oil of the latter is usually extracted with solvents such as petroleum ether or n-hexane in order to obtain rose concrete. Rose concrete, the result of solvent extraction, is mainly composed of fragrance-related substances, and contains large quantities of paraffins, fatty acids, fatty acid methyl esters, di- and tri-terpenic

therefore processed through Solvent-Free Microwave Extraction (SFME) as upcoming extraction technique for the extraction of essential oil to eliminate non-volatile compounds such as paraffins, shorten extraction time, reduce organic solvent consumption, improve extraction yield, enhance quality of the extract, prevent pollution and reduce sample preparation costs [1].

[1] M. Mohamadi, T. Shamspur, A. Mostafavi, Journal of Essential Oil Research 2012, 25, 55-61.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from Rosa Damascena Mill. can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-Free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHD techniques. (See Cookbook for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Fresh Rosa Damascena Mill., SFME						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1kW	2.1kW		
Small	400	400	•		0.4	0.1
Medium	1264	1264		•	1.8	0.14
Large	2978	1800		•	3.9	0.13

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

## Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows you to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

In this study, we propose state-of-the-art processes for extraction of essential oils from Dry Rosa

Damascena Mill. through SFME and MWHD. It is the unique modern concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

Compound	RI <sup>a</sup>	SFME (% <sup>b</sup> )
<b>Monoterpenes</b>		<b>0.6</b>
$\alpha$ -Pinene	858	0.2
Sabinene	896	tr <sup>c</sup>
$\beta$ -Pinene	898	0.1
$\beta$ -Myrcene	916	0.1
$\alpha$ -Terpinene	937	tr
Limonene	949	tr
$\gamma$ -Terpinene	977	tr
$\alpha$ -Terpinolene	1003	tr
<b>Oxygenated monoterpenes</b>		<b>49.3</b>
Rosefuran	1009	tr
Linalool	1017	0.4
<i>trans</i> -Rose oxide	1023	0.3
Citronellal	1056	tr
Neroloxide	1061	tr
4-Terpineol	1081	0.1
$\beta$ -Citronellol	1157	34.7
Geraniol	1173	9.1
Eugenol	1262	1.5
Geranyl acetate	1596	0.6
Methyl eugenol	1309	2.3
Neryl acetate	1926	0.2
<b>Sesquiterpenes</b>		<b>4.1</b>
$\beta$ -Bourbonene	1302	0.2
<i>trans</i> -caryophyllene	1326	0.5

Table 2. Qualitative and quantitative composition of rose essential oils obtained by SFME

Compound	RI <sup>a</sup>	SFME (% <sup>b</sup> )
$\alpha$ -Guaiene	1340	0.4
$\alpha$ -Humulene	1348	0.5
Germacrene D	1368	1.4
$\beta$ -Selinene	1369	0.1
Caryophyllen (1I)	1378	0.1
$\alpha$ -Selinene	1379	-
$\delta$ -Guaiene	1381	0.4
( <i>E,E</i> )- $\alpha$ -Farnesene	1385	0.1
$\delta$ -Cadinene	1394	tr
Ledene	1518	0.3
<b>Oxygenated sesquiterpenes</b>		<b>2.2</b>
Elmol	1419	tr
Nerolidol	1437	0.1
Caryophyllene oxide	1446	tr
$\gamma$ -Eudesmol	1505	0.1
$\beta$ -Eudesmol	1516	0.2
<i>cis</i> -Farnesol	1555	1.8
<b>Hydrocarbons</b>		<b>31.4</b>
Pentadecane	1388	0.2
Hexadecane	1492	0.1
8-Heptadecene	1547	0.2
Heptadecane	1569	1.8
<i>cis</i> -9-Tricosene	1630	0.2
Octadecane	1669	0.4
( <i>Z</i> )-5-Nonadecene	1737	3.5
Nonadecane	1765	15.1
( <i>E</i> )-9-Eicosene	1811	0.6
Eicosane	1836	1.2
9-Nonadecene	1882	0.2
1-Nonadecene	1894	0.2
<i>n</i> -Heneicosan	1910	6.0
Docosan	1972	0.1
<i>cis</i> -9-Tricosene	2029	0.2
Tricosane	2039	0.9
Tetracosane	2099	tr
Pentacosane	2160	0.2
<i>n</i> -Heptacosane	2274	0.2
<b>Other oxygenated compounds</b>		<b>0.6</b>
Phenylethyl alcohol	813	0.4
Heptanal	1011	tr

Table 2 (continued).

Compound	RI <sup>a</sup>	SFME (% <sup>b</sup> )
Tetradecanal	1024	0.1
Nonanal	1485	tr
<b>Total oxygenated compounds</b>		<b>52.1</b>
<b>Total non-oxygenated compounds</b>		<b>36.1</b>
<b>Extraction time (min)</b>		<b>35</b>

<sup>a</sup>Retention indices relative to C6–C27 *n*-alkanes on HP-1MS column.

<sup>b</sup>%, relative percentage obtained on HP-1 column using GC/MS detector.

<sup>c</sup>tr, 60.05. \*Significant at  $p < 0.05$  based on *F*-value determined by analysis of variance (ANOVA). \*\*Significant at  $p < 0.01$  based on *F*-value determined by ANOVA. CV, coefficient of variation.

Table 2 (continued).



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX15 - ROSEMARY FLAVORS

---

Rosemary Microwave  
Hydrodiffusion and Gravity (MHG)



### Introduction

Rosemary (*Rosmarinus officinalis* L.) is a perennial herb with fragrant evergreen needle-like leaves. It is native to the Mediterranean region and it has been cultivated for a long time. It belongs to the Lamaiaceae family, which comprises up to 200 genera and about 3500 species. The leaves are evergreen, with dense short woolly hairs. Rosemary has been a significant herb since antiquity, although rosemary is more familiar to contemporary Westerners as a kitchen herb used to add a spicy or slightly medicinal flavor to some foods, it was traditionally used as an antiseptic, astringent, and food preservative before the invention of refrigeration. Rosemary's antioxidant properties are still used to extend the shelf life of prepared foods. Rosemary is also known medicinally for its powerful antioxidant activity, its antibacterial and antimutagenic properties, and as a chemopreventive agent. Besides the therapeutical application, the essential oil is widely applied in the cosmetic industry producing various Cologne waters, bathing essences, hair lotions and shampoos and as a component of disinfectants and insecticides [1].

### Why to choose Microwave Flavor set-up?

The patented and revolutionary Microwave Hydrodiffusion and Gravity (MHG) system paves the way to new flavoring products which were impossible to be obtained with the ancient extraction concepts. Percolation, solvent extraction etc... were inefficient and environmental-unfriendly methods of flavor extraction. MHG is going to improve the efficiency and the quality of flavoring products.

- New natural flavors
- Fast extraction
- No thermal degradation
- No solvent

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance were obtained through the MHG techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and experimental procedure

The MHG technique is suitable for both dry and fresh raw material, see the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Fresh Rosemary (MHG)						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Total flavour extract [mL]	Total flavor extract yield [%]
			1kW	2.1kW		
Small	500	500	•		267	53.3
Medium	1580	1580		•	855	54.1
Large	3720	1800		•	2001	53.8

Dry Rosemary (MHG)						
Reactor	Weighted dry soaked material [g]	Power [W]	Chiller		Total flavour extract [mL]	Total flavor extract yield [%]
			1 kW	2.1kW		
Small	500	500	•		180	36
Medium	1580	1580		•	540	34.2
Large	3720	1800		•	1330	35.8

\*Time, Power

≤ 1800 g: Power(W) = Weight(g) for 20 min.

> 1800g: Power = 1800W for 40min

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1



**Important remarks**

Please take into account that the interaction between microwaves and raw material is based on the water content. The optimized method (power and time) depends on the type of Rosemary, mainly to the water content of Rosemary itself. Please use the reported parameter as general application note to start to optimize your own method. Be careful, that using an excess power might cause burning of your sample.

Please take care to seal properly the glass reactor during the installation of the flavors set-up according to the manual, to avoid loss of vapor during extraction.

**Conclusion**

A newly and cleaner design process for extraction of flavors was developed in this study: MHG. This green process has been studied and tested using Rosemary. This new system was developed to date indicate that microwave extraction process of flavors offers important advantages over antiquated extraction techniques, namely, shorter extraction times, less energy consuming, lower costs as well as new flavoring products. The MHG system offers furthermore the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX16 - ROSEMARY FRAGRANCES

---

Rosemary Solvent-Free Microwave Extraction (SFME) and Microwave Hydrodistillation (MWHD)



### Introduction

Rosemary (*Rosmarinus officinalis* L.) is an important herb on the world food and aromatherapy market. Natural antioxidants such as those present in rosemary essential oil may be an alternative source for compounds capable of protecting lipids in foods. The essential oil secreted by glandular trichomes is mainly located in leaves. Rosemary essential oil is also used as an antibacterial and antifungal agent. Nevertheless, it has been noted that these activities often depend on the origin of the rosemary plant and the method of extraction. Since both of these quality parameters can greatly influence the chemical composition of rosemary oil, Solvent-Free Microwave Extraction (SFME) as upcoming extraction techniques have been reported for the extraction of fragrances from Rosemary [2].

[2] N. Tigrine-Kordjani, B. Meklati, F. Chemat, International Journal of Aromatherapy 2006, 16, 141-147.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from Rosemary can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-Free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHD techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 2*).

Fresh Rosemary, SFME						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1kW	2.1kW		
Small	500	500	•		3.2	0.63
Medium	1580	1580		•	10.3	0.65
Large	3720	1800		•	26.8	0.72

Dry Rosemary, MWHD						
Reactor	Weighted dry soaked material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1 kW	2.1kW		
Small	500	500	•		1.9	0.38
Medium	1580	1580		•	6.0	0.38
Large	3720	1800		•	14.9	0.4

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 2

## Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows you to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

In this study, we propose state-of-the-art processes for the extraction of essential oils from rosemary through SFME and MWHD. It is the unique modern concept of

the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.

N°	Compound <sup>a</sup>	SFME (%)	R <sup>b</sup>	R <sup>c</sup>
	Monoterpene hydrocarbons	68.56		
1	Tricyclene	0.26	921	1011
2	$\alpha$ -Pinene	44.05	936	1023
3	Camphene	6.14	951	1103
4	Verbenene	0.77	955	1121
5	$\beta$ -Pinene	2.61	980	1109
6	Myrcene	1.94	995	1149
7	$\alpha$ -Phellandrene	0.31	995	1165
8	$\gamma$ -3-Carene	0.08	1014	1290
9	$\alpha$ -Terpinene	0.86	1020	1083
10	para-Cymene	1.27	1025	1250
11	Limonene	5.48	1030	1206
12	$\gamma$ -Terpinene	3.08	1052	1251
13	Terpinolene	1.71	1092	1287
	Oxygenated monoterpenes	24.87		
14	Linalool	2.00	1106	1538
15	$\alpha$ -Campholenal	1.24	1122	1471
16	Camphor	7.82	1149	1514
17	Pinocarvone	1.33	1160	1548
18	Borneol	2.57	1173	1679
19	Terpin-4-ol	2.07	1184	1590
20	$\alpha$ -Terpineol	0.77	1198	1677
21	Verbenone	6.37	1207	1696
22	Geraniol	0.70	1279	1828

Table 3. Chemical composition of *Rosmarinus officinalis* essential oils obtained by SFME.

N°	Compound <sup>a</sup>	SFME (%)	RI <sup>b</sup>	RI <sup>c</sup>
	Sesquiterpene hydrocarbons	1.91		
23	E-Caryophyllene	0.95	1425	1470
24	$\alpha$ -Humulene	0.42	1450	1657
25	$\gamma$ -Curcumene	0.04	1469	1738
26	$\beta$ -Bisabolene	0.43	1508	1714
27	$\beta$ -Sesquiphellandrene	0.07	1519	1776
	Oxygenated sesquiterpenes	0.26		
28	Caryophyllene oxide	0.10	1570	1977
29	$\alpha$ -Bisabolol	0.16	1684	2022
	Other oxygenated compounds	1.03		
30	Bornyl acetate	0.81	1263	1579
31	Methyl eugenol	0.12	1397	2032
32	Z-Methyl jasmonate	0.10	1635	2349
	Extraction time (min)	45		
	Total oxygenated compounds	26.16		
	Total non-oxygenated compounds	70.47		

<sup>a</sup> Essential oil compounds sorted by chemical families and percentages calculated by GC–FID on non-polar HP5MSTM capillary column.

<sup>b</sup> Retention indices relative to C<sub>5</sub>–C<sub>28</sub> n-alkanes calculated on non-polar HP5MS™ capillary column.

<sup>c</sup> Retention indices relative to C<sub>5</sub>–C<sub>28</sub> n-alkanes calculated on polar Carbowax™-PEG capillary column.

*Table 3 (continued).*



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX17 - THYMUS FRAGRANCES

---

Thymus Vulgaris Solvent-Free  
Microwave Extraction (SFME) and  
Microwave Hydrodistillation (MWHD)



### Introduction

Thymus Vulgaris L. (common thyme), a member of the Labiatae family, is an aromatic/medicinal plant of increasing economic importance for North America, Europe, North Africa and Asia. Thyme is one of many aromatic plants that has been utilized in a variety of food products to provide a flavor specific to this herb. Studies indicating the antiseptic, carminative, antimicrobial, and antioxidative properties of thyme have also been published. From the medicinal point of view, thyme has been used as a culinary herb and also as a herbal medicine. The essential oils of thyme are responsible for the typical spicy aroma of the plant leaves. These oils are stored in glandular peltate trichomes situated on both sides of the leaves. Results published on the chemical composition of thyme oil revealed that most of the oil was produced from flowering plants. In the plant's life cycle, the oil production is usually at its highest level during this period. For *T. vulgaris* and *T. pulegioides*, such finding was reported in the early 1960s <sup>[1]</sup>.

### Why Choose Microwave Fragrances Set-up?

The standard method is the Clavenger method, which was published for the first time in 1928. According to that method, the essential oil from *Thymus Vulgaris* can be extracted by hydrodistillation or steam distillation. These techniques take several hours of heating which may cause degradation of thermolabile compounds present in the starting plant material and therefore odor deterioration. The patented and innovative Microwave Hydrodistillation (MWHD) and Solvent-Free Microwave Extraction (SFME) techniques allow the production of essential oils with higher quality.

- High quality fragrances
- No thermal degradation
- Fast extraction
- High purity, no artifacts

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality fragrance was obtained through the SFME or MWHD techniques (see the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The SFME and the MWHD techniques are respectively suitable for fresh and dry raw materials. See the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Fresh <i>Thymus Vulgaris</i> , SFME						
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1kW	2.1kW		
Small	400	400	•		1.3	0.32
Medium	1264	1264		•	4.4	0.35
Large	2978	1800		•	11	0.37

Dry <i>Thymus Vulgaris</i> , MWHD						
Reactor	Weighted dry soaked material [g]	Power [W]	Chiller		Volatile fraction [mL]	Yield [%]
			1 kW	2.1kW		
Small	400	400	•		0.6	0.16
Medium	1264	1264		•	1.9	0.15
Large	2978	1800		•	5.4	0.18

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g).

> 1800g: Power = 1800W

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

## Important Remarks

The system is developed with an automatic recirculation of the distilled water. This allows you to manage extraction power and time to match your own specific requirements. Please make sure to properly seal the glass reactor during the installation of the fragrances set-up according to the manual, to avoid loss of vapor during extraction.

## Conclusion

In this study, we propose state-of-the-art processes for the extraction of essential oils from *Thymus Vulgaris* through SFME and MWHD. It is the unique modern

concept of the antiquated Clavenger method, highly accelerating the isolation process, without causing changes in the volatile oil composition. The efficiency of the new techniques SFME and MWHD are considerably higher than the conventional procedure, if we take into account short distillation times required, cost and energy used and cleanliness of the process. An added-value feature is the possibility to work with different scalar matrices amounts due to three different volumes of the reactor vessels (small, medium, large) complying with a high range of extraction-scale needs.

N°	RT <sup>b</sup> (min)	Compound	I <sub>k</sub> <sup>c</sup>	Relative peak area <sup>a</sup> (%) SFME
1	16.05	$\alpha$ -Thujene	930	0.53 $\pm$ 0.02
2	16.53	$\alpha$ -Pinene	938	0.86 $\pm$ 0.01
3	17.31	Camphene	952	0.53 $\pm$ 0.01
4	19.04	1-Octen-3-ol	983	2.64 $\pm$ 0.31
5	19.60	$\beta$ -Myrcene	993	1.30 $\pm$ 0.17
6	19.91	3-Octanol	999	0.19 $\pm$ 0.03
7	20.39	$\alpha$ -Phellandrene	1008	0.18 $\pm$ 0.01
8	20.71	$\Delta$ -3-Carene	1014	0.09 $\pm$ 0.01
9	21.12	$\alpha$ -Terpinene	1021	1.73 $\pm$ 0.14
10	21.73	<i>p</i> -Cymene	1033	17.57 $\pm$ 0.78
11	21.99	1,8-Cineole	1037	1.31 $\pm$ 0.12
12	23.55	$\gamma$ -Terpinene	1066	8.54 $\pm$ 0.02
13	23.91	Trans-Sabinene hydrate	1073	0.94 $\pm$ 0.05
14	24.95	Terpinolene	1093	0.27 $\pm$ 0.05
15	25.59	Linalool	1105	2.43 $\pm$ 0.27
16	29.20	Borneol	1176	1.11 $\pm$ 0.21
17	29.69	Endo-Borneol	1185	1.41 $\pm$ 0.21
18	29.80	Terpinen-4-ol	1187	0.63 $\pm$ 0.16
19	31.20	$\alpha$ -Terpineol	1216	0.17 $\pm$ 0.00
20	32.34	Methyl thymylether	1240	0.14 $\pm$ 0.07
21	34.32	Geraniol	1281	0.39 $\pm$ 0.07
22	35.90	Thymol	1315	40.20 $\pm$ 3.03
23	36.50	Carvacrol	1328	6.84 $\pm$ 0.68

Table 2. Chemical compositions of essential oils obtained by solvent-free microwave extraction (SFME) of thyme aerial parts.



N°	RT <sup>b</sup> (min)	Compound	I <sub>K</sub> <sup>c</sup>	Relative peak area <sup>a</sup> (%) HD
24	38.22	Thymyl acetate	1366	0.16 ± 0.03
25	41.36	β-Caryophyllene	1438	2.86 ± 0.27
26	42.72	α-Humulene	1470	0.64 ± 0.22
27	42.98	Geranyl acetate	1477	0.42 ± 0.09
28	45.44	Δ-Cadinene	1536	0.40 ± 0.05
29	48.17	Caryophyllene oxide	1604	1.42 ± 0.21
Total peak area (%)			94.31	95.91

<sup>a</sup> Mean ± SD (n = 2).

<sup>b</sup> Retention time.

<sup>c</sup> Kovats Retention Index (I<sub>K</sub>) relative to C<sub>9</sub>-C<sub>18</sub> n-alkanes on the HP-5MS column.

*Table 2 (continued).*



**MILESTONE**  
H E L P I N G  
C H E M I S T S

# APPLICATION REPORT

## EX18 - STRAWBERRY FLAVORS

---

### Strawberry Microwave Hydrodiffusion and Gravity (MHG)



#### **Introduction**

Extraction, notably from fruits for production of juices or specific ingredients, is a highly studied subject. Microwave Hydrodiffusion (MWH) can be used for the extraction of flavors from fruits or for the production of various extracts. Fruits and vegetables are known for their health beneficial properties, however the nutritional recommendations have proven hard to reach by the general population, among others due to lack of convenience. Fruit juice is viewed by the general public as a convenient alternative to fresh fruits and the fruit juice production presents a high economical attraction (Étiévant et al., 2010). The health beneficial impact of fruits and vegetables are connected to their content in micronutrients and antioxidants, so the density of fruit and vegetable products in micronutrients and antioxidants has become a preoccupation of the food industry. A challenge is therefore to produce fruit juices with high nutritional densities, for example by using fruits rich in antioxidants.

Aurélie Cendres, F. C. (2012).  
Comparison between microwave hydrodiffusion and pressing for plum juice. *LWT - Food Science and Technology*, 229-237.

### Why Choose Microwave Flavor Set-up?

The patented and revolutionary Microwave Hydrodiffusion and Gravity (MHG) system paves the way to new flavoring products which were impossible to be obtained with the ancient extraction concepts. Percolation, solvent extraction, etc. were inefficient and not environmentally-friendly methods of flavor extraction. MHG improves the efficiency and the quality of flavoring products.

- New natural flavors
- Fast extraction
- No thermal degradation
- No solvent

### Instrumentation and Principles of Operation

A very efficient extraction process can be achieved thanks to the selective heating of microwaves to materials through molecular interactions with the electromagnetic field via conversions of electromagnetic energy into thermal energy. The high quality flavor was obtained through MHG techniques. (See the "Microwave Extraction Techniques" section for theory and principle).

### Results and Experimental Procedure

The MHG technique is suitable for both dry and fresh raw material, see the "Quick start guide" for a list of easy and sequential setting-up operations (*Table 1*).

Fresh Strawberry							
Reactor	Weighted fresh raw material [g]	Power [W]	Chiller		Total flavor extract [mL]	Volatile fraction [mL]	Total flavor extract yield [%]
			1kW	2.1kW			
Small	610	610	•		350	0	57.2
Medium	1920	1800		•	1100	0	57.0
Large	4540	1800		•	2600	0	57.1

#### Time, Power

The extractions were carried out till complete recovery of the fragrance

≤ 1800 g: Power(W) = Weight(g) for 20 min.

> 1800g: Power = 1800W for 40min

Chiller settings:

≤ 900g, 1 kW Chiller

> 900g, 2.1 kW Chiller

Table 1

### **Important Remarks**

Please take into account that the interaction between microwaves and raw material is based on the water content. The optimized method (power and time) depends on the type of strawberries, mainly to the water content of the strawberry itself. Please use the reported parameter as general application note to start to optimize your own method. Be careful, that using an excess power might cause burning of your sample. Please make sure to properly seal the glass reactor during the installation of the flavors set-up according to the manual, to avoid loss of vapor during extraction.

### **Conclusion**

A new and cleaner design process for extraction of flavors was developed in this study: MHG. This green process has been studied and tested using strawberries. This new system was developed to indicate that the microwave extraction process of flavors offers important advantages over antiquated extraction techniques, namely, shorter extraction times, less energy consuming, lower costs as well as new flavoring products. Furthermore, the MHG system offers the possibility to work with scalar amounts of sample due to three different reactor vessels (small, medium, large), complying with a high range of extraction-scale needs.